



1289 South Rd, St Mary's 5042 South Australia
 Telephone: 1300 266 730
 Fax: 08 8177 2145
 Email: Info@sandgoanna.com.au
 Website: www.sandgoanna.com.au

Sand Goanna
 Quality Industry Products

ABN: 46 994 171 449

Automatic External Defibrillator (Trainer AED)



2 Colours

Call Today:

1300 266 730



Special...
\$270.00

10 preprogrammed scenarios in The AED Trainer:

Scenarios	Description	Operation
0	Ventricular Fibrillation- With single-shock conversion	Shockable rhythm 1 shock Non-shockable rhythm
1	Ventricular Fibrillation- With multiple shocks required for conversion	Shockable rhythm 4 shocks Non-shockable rhythm
2	Troubleshooting Skills- Defibrillation pads	Poor pad connection Shockable rhythm 1 shock Non-shockable rhythm
3	Ventricular Fibrillation- With two shocks required for conversion	Shockable rhythm 2 shock Non-shockable rhythm
4	Non-Shockable Rhythm	Non-shockable rhythm throughout
5	Ventricular Fibrillation- With two shocks required for conversion	Shockable rhythm 2 shocks Non-shockable rhythm
6	Ventricular Fibrillation- With three shocks required for conversion	Shockable rhythm 3 shocks Non-shockable rhythm
7	Troubleshooting Skills- Defibrillation pads with two shocks required for conversion	Poor pad connection Shockable rhythm 2 shocks Non-shockable rhythm
8	Ventricular Fibrillation	Ventricular fibrillation throughout
9	Troubleshooting Skills- Motion, low battery, with one shock required for conversion	Motion artifacts Shockable rhythm

This Trainer is compatible with the new guidelines from AHA2009 scenarios & operations.

The AED trainer is designed to prepare emergency responders to use the real clinical AED.

This trainer provides a variety of simulations, or training scripts, to help responders become familiar with the AED and allow them to demonstrate the basic skills necessary to use the AED in an emergency.

Features

AED trainer provides simulated shock delivery. It has no high-voltage capabilities, ensuring safety during training.

Note: Remember handling safety when using working defibrillators that give a high-voltage shock.

The AED trainer is designed for use with human body. This link technology of XFT company can give responses to pad placement on human body. It can also be used with any other manikin.

The AED Trainer is pre-configured with 10 training scenarios that simulate realistic sudden cardiac arrest episodes. See the scenario descriptions in page 8 of manual. These scenarios are compatible with training programs developed by nationally and internationally recognized responder programs. Check with your local medical director and Emergency Medical System for additional training requirements, protocols, and materials.